



## **Optimizing Performance of Transcreener Fluorescence Polarization Assays with the Feyond-A500 Multi-Mode Microplate Reader**

Transcreener technology is a universal, high-throughput biochemical assay platform based on nucleotide detection.


The assay is based on the detection of nucleotide diphosphates (ADP / GDP).

Nucleotide diphosphates are formed by thousands of kinases, many of which catalyze covalently regulated reactions.

These reactions are central to cellular signaling and are of great value in drug discovery.

## Principle

The reagents for all the assays are a far-red tracer bound to a highly specific monoclonal / polyclonal antibody. An enzymatic reaction generates diphosphates or monophosphates, which displace the tracer from the antibody-quencher conjugate. This results in the generation of a signal due to an increase in rotational freedom of the tracer, detected as a decrease in polarization.



## Verification Standards

Prepare a 10  $\mu\text{M}$  ATP / ADP standard curve to simulate the enzyme reaction. Starting with 10  $\mu\text{M}$  ATP, increase the amount of ADP added and decrease ATP proportionally, maintaining the total adenine nucleotide concentration at 10  $\mu\text{M}$ . At a 10% conversion rate of 10  $\mu\text{M}$  ATP,  $Z' > 0.7$  and  $\Delta\text{mP} > 120$ .

## Materials and Methods

- Feyond-A500 Multi-Mode Microplate Reader (Hangzhou Allsheng)
- Transcreener® ADP<sup>2</sup> FP Assay (Code: No.3010-1K)
- ATP / ADP Mixture In Buffer (Constant Adenine Concentration: 10  $\mu\text{M}$ )
- Corning 96-Well Black Polystyrene Plate (Code: No. 3915)

Concentration (%)	ATP	ADP
100	0	100
75	25	75
50	50	50
25	75	25
15	85	15
10	90	10
7.5	92.5	7.5
5	95	5
3	97	3
2	98	2
1	99	1
0	100	0

Table 1 ADP / ATP Standard Curve Preparation (10  $\mu\text{M}$ )

Parameter	Fluorescence Polarization
EX	624-40 nm
EM	692-40 nm
G-factor	1.06
Number	150
PMT gain	Auto
Integration time	40 $\mu\text{s}$

Table 2 Instrument Parameter Settings





Figure 1 FP Measurement Mode Selection




Figure 2 Instrument Parameter Settings

## Results

As the ratio of ADP to ATP increases, the proportion of bound tracer vs. free tracer decreases, resulting in an overall decrease in mP values.



A: Z' and  $\Delta mP$  values observed in a standard curve mimic the conversion of 10  $\mu M$  ATP to ADP.



B: Zoomed view of the 0-2.5  $\mu M$  ADP section of the standard curve shows the Z' validation minimal qualification data (blue dashed line) and  $\Delta mP$  validation minimal qualification data (red dashed line). The 10% ATP conversion validation point is also indicated (vertical black dotted line).

## Conclusion

The Feyond-A500 Multi-Mode Microplate Reader passed the validation criteria for the Transcreener ADP<sup>2</sup> FP assay. The filter-based measurement results showed a Z' value of 0.79 (standard: Z' value > 0.70 at 10% conversion at 10  $\mu M$  ATP) and a  $\Delta mP$  of 140 (standard:  $\Delta mP$  > 120 mP at 10% conversion at 10  $\mu M$  ATP).

## HANGZHOU ALLSHENG INSTRUMENTS CO., LTD.

Building 9 No.7 of Zhuantang Science and Technology Economic Zone,  
Xihu District, Hangzhou City, 310024 Zhejiang, P.R. China

Tel: +86-571-88859758

Fax: +86-571-87205673

✉ [info@allsheng.com](mailto:info@allsheng.com)

🌐 [www.allsheng.com](http://www.allsheng.com)

Version Date: 2025.09

